

# TMSS10ListContainer

TMSS10ListContainer lists containers and any assigned Volumes

## Synopsis

```
TMSS10ListContainer [-A attribute index] [-S logon string] [-V volume filter]
[-c configuration file] [-s sync friendly output]
```

## Options

- -A Extended Attribute index (default = 0).
- -S The Logon string argument sets the [Server Logon Information](#).
- -V [Volume Filter](#) which can be used to filter the traversal of Customer, Media, Volume, repository and flag values.
- -c The config file argument specifies the path to the Batch Configuration File. The Batch Configuration File allows you to perform advanced traversal of [Customer](#), [Media](#) and [Volumes](#), perform advanced Barcode pre and post processing and route program output to syslogd.
- -s Produces output in a sync friendly format, including listing container-id relisted for each Volume entry.

If an argument value starts with the value FI:, the value for that attribute will be read from the file name specified immediately after it. For example, if a file named ACME\_volumes has the content ACME\_\*.\*, specifying FI:/etc/ACME\_volumes will have the same effect as -V ACME.\*.\*. Using FI: to refer to a file for an attribute value works for all attributes.

If the pattern being matched starts with \*LIST: then rather than matching the literal value, each record in the file specified as \*LIST:filename will be tested. If one matches the match will succeed.

2025/01/21 22:08

[include page](#)

To include multiple Attributes in the report, add the required Attribute indexes separated by a comma (eg. -A 1,2,5).

2025/01/21 22:08

[include page](#)

## Technical Support

The TapeTrack Software is commercially supported by a full-time help desk staff.

If you are experiencing problems or want some advice on how to configure or use the product please see the [Accessing Technical Support](#) page.

Files

**stdout:** Output Reports.  
**stderr:** Diagnostic messages.

Example

Standard format output

Command line syntax

```
TMSS10ListContainer -S user:-password@server
```

Output

|                                                                               |                   |               |      |                     |
|-------------------------------------------------------------------------------|-------------------|---------------|------|---------------------|
| System: TapeTrack                                                             |                   |               |      |                     |
| Report ID: 123                                                                |                   |               |      |                     |
| Report Description: Container List                                            |                   |               |      |                     |
| Production Time: Friday, June 23, 2023 - 14:35:23 (AUS Eastern Standard Time) |                   |               |      |                     |
| Report Width: 79 bytes                                                        |                   |               |      |                     |
| Seq.                                                                          | Barcode           | Capacity Used | Item | Content             |
| -----                                                                         |                   |               |      |                     |
| --                                                                            |                   |               |      |                     |
| 1                                                                             | US01.CONT.CONT001 | 8             |      | 1 US01.LT0.000001L6 |
|                                                                               |                   |               |      | 2 US01.LT0.000002L6 |
|                                                                               |                   |               |      | 3 US01.LT0.000003L6 |
|                                                                               |                   |               |      | 4 US01.LT0.000004L6 |
|                                                                               |                   |               |      | 5 US01.LT0.000005L6 |
|                                                                               |                   |               |      | 6 US01.LT0.000006L6 |
|                                                                               |                   |               |      | 7 US01.LT0.000007L6 |
|                                                                               |                   |               |      | 8 US01.LT0.000008L6 |
| 2                                                                             | US01.CONT.CONT002 | 7             |      | 1 US01.LT0.000011L6 |
| 3                                                                             | US01.CONT.CONT003 |               |      | 2 US01.LT0.000012L6 |
|                                                                               |                   |               |      | 3 US01.LT0.000013L6 |
|                                                                               |                   |               |      | 4 US01.LT0.000014L6 |
|                                                                               |                   |               |      | 5 US01.LT0.000015L6 |
|                                                                               |                   |               |      | 6 US01.LT0.000018L6 |
|                                                                               |                   |               |      | 7 US01.LT0.000021L6 |
| 4                                                                             | US01.CONT.CONT004 | 5             |      | 1 US01.LT0.000025L6 |
|                                                                               |                   |               |      | 2 US01.LT0.000027L6 |
|                                                                               |                   |               |      | 3 US01.LT0.000030L6 |
|                                                                               |                   |               |      | 4 US01.LT0.000031L6 |
|                                                                               |                   |               |      | 5 US01.LT0.000033L6 |

```
5 US01.CONT.CONT005
6 US01.CONT.CONT006
7 US01.CONT.CONT007
8 US01.CONT.CONT008
9 US01.CONT.CONT009
10 US01.CONT.CONT010
-----
--
```

Sync friendly format output

Command line syntax

```
TMSS10ListContainer -S user:-password@server
```

Output

System: TapeTrack  
Report ID: 123  
Report Description: Container List  
Production Time: Friday, June 23, 2023 - 14:37:25 (AUS Eastern Standard Time)  
Report Width: 79 bytes

| Seq. | Barcode             | Capacity Used | Item | Content           |
|------|---------------------|---------------|------|-------------------|
| --   |                     |               |      |                   |
|      |                     |               |      |                   |
|      | 1 US01.CONT.CONT001 | 8             | 1    | US01.LT0.000001L6 |
|      | 1 US01.CONT.CONT001 | 8             | 2    | US01.LT0.000002L6 |
|      | 1 US01.CONT.CONT001 | 8             | 3    | US01.LT0.000003L6 |
|      | 1 US01.CONT.CONT001 | 8             | 4    | US01.LT0.000004L6 |
|      | 1 US01.CONT.CONT001 | 8             | 5    | US01.LT0.000005L6 |
|      | 1 US01.CONT.CONT001 | 8             | 6    | US01.LT0.000006L6 |
|      | 1 US01.CONT.CONT001 | 8             | 7    | US01.LT0.000007L6 |
|      | 1 US01.CONT.CONT001 | 8             | 8    | US01.LT0.000008L6 |
|      | 2 US01.CONT.CONT002 |               |      |                   |
|      | 3 US01.CONT.CONT003 | 7             | 1    | US01.LT0.000011L6 |
|      | 3 US01.CONT.CONT003 | 7             | 2    | US01.LT0.000012L6 |
|      | 3 US01.CONT.CONT003 | 7             | 3    | US01.LT0.000013L6 |
|      | 3 US01.CONT.CONT003 | 7             | 4    | US01.LT0.000014L6 |
|      | 3 US01.CONT.CONT003 | 7             | 5    | US01.LT0.000015L6 |
|      | 3 US01.CONT.CONT003 | 7             | 6    | US01.LT0.000018L6 |
|      | 3 US01.CONT.CONT003 | 7             | 7    | US01.LT0.000021L6 |
|      | 4 US01.CONT.CONT004 | 5             | 1    | US01.LT0.000025L6 |
|      | 4 US01.CONT.CONT004 | 5             | 2    | US01.LT0.000027L6 |
|      | 4 US01.CONT.CONT004 | 5             | 3    | US01.LT0.000030L6 |
|      | 4 US01.CONT.CONT004 | 5             | 4    | US01.LT0.000031L6 |
|      | 4 US01.CONT.CONT004 | 5             | 5    | US01.LT0.000033L6 |

```
5 US01.CONT.CONT005
6 US01.CONT.CONT006
7 US01.CONT.CONT007
8 US01.CONT.CONT008
9 US01.CONT.CONT009
10 US01.CONT.CONT010
```

```
-----
--
```

## Changelog

Cannot load rss feed.

From:

<https://docs.tapetrack.com/> - **TapeTrack Documentation**

Permanent link:

<https://docs.tapetrack.com/cli/tmss10listcontainer?rev=1687495132>

Last update: **2025/01/21 22:07**

