

Media Options : Relaxed Barcode Matching

The [Media](#) Option Relaxed Barcode Matching will, when a scanned barcode is not found, attempt to match the Barcode with any [Volumes](#) with a L|T|J suffix.

Example

If a barcode 100123 is scanned and not found, Relaxed Barcode Matching will then attempt to match against the barcodes L|T|J suffix, such as 100123L4, 100234T1 or 100234J2.

Options

- True: Enables Relaxed Barcode Matching.
- False: Disables Relaxed Barcode Matching.

[options](#), [media options](#), [barcode](#)

From:

<https://docs.tapetrack.com/> - TapeTrack Documentation

Permanent link:

https://docs.tapetrack.com/options/media/relaxed_barcode_matching?rev=1557883618

Last update: 2025/01/21 22:07

