


# Relaxed Barcode Matching

The [Media Option](#) Relaxed Barcode Matching will, when a scanned Barcode is not found, attempt to match the Barcode with any [Volumes](#) with a L|T|J suffix.



Edit Media: LTO Cartridge (6/0)

Identity Volume-ID Building Options

|                           |         |
|---------------------------|---------|
| Hide Media Type           | False   |
| Internal                  | False   |
| Maximum Volumes           | 0       |
| Maximum Repositories      | 0       |
| Allowable Volume-IDs      |         |
| Default Media Type        | False   |
| Auto Simple Management    | False   |
| Auto Container            | False   |
| Stock Controlled          | False   |
| Auto Stock Selection      | False   |
| Relaxed Barcode Matching  | False   |
| No Move-back Warning      | True    |
| Lock Volumes              | False   |
| No HealthCheck            | False   |
| Show BuiltIn Repositories | False   |
| Remove barcode suffix     |         |
| Barcodes use Modulo43     | False   |
| Default Repository        | Library |

**Relaxed Barcode Matching**  
Match scanned barcode against barcode with L|T|J suffix when not found

Save Cancel Help

## Options

- True: Enables Relaxed Barcode Matching.
- False: Disables Relaxed Barcode Matching.

# Examples

If a barcode 100123 is scanned and not found, Relaxed Barcode Matching will then attempt to match against the barcodes L|T|J suffix, such as 100123L4, 100234T1 or 100234J2.

[options](#), [media options](#), [barcode](#)

From:  
<https://docs.tapetrack.com/> - **TapeTrack Documentation**

Permanent link:  
[https://docs.tapetrack.com/options/media/relaxed\\_barcode\\_matching?rev=1573428224](https://docs.tapetrack.com/options/media/relaxed_barcode_matching?rev=1573428224)

Last update: **2025/01/21 22:07**

