


Relaxed Barcode Matching

The [Media Option](#) Relaxed Barcode Matching will, when a scanned Barcode is not found, attempt to match the Barcode with any [Volumes](#) with a L|T|J suffix.



Edit Media: LTO Cartridge (6/0)

Identity Volume-ID Building Options

Hide Media Type	False
Internal	False
Maximum Volumes	0
Maximum Repositories	0
Allowable Volume-IDs	
Default Media Type	False
Auto Simple Management	False
Auto Container	False
Stock Controlled	False
Auto Stock Selection	False
Relaxed Barcode Matching	False
No Move-back Warning	True
Lock Volumes	False
No HealthCheck	False
Show BuiltIn Repositories	False
Remove barcode suffix	
Barcodes use Modulo43	False
Default Repository	Library

Relaxed Barcode Matching
Match scanned barcode against barcode with L|T|J suffix when not found

Save Cancel Help

Options

- True: Enables Relaxed Barcode Matching.
- False: Disables Relaxed Barcode Matching.

Examples

If a Barcode 100123 is scanned and not found, Relaxed Barcode Matching will then attempt to match against the Barcodes L|T|J suffix, such as 100123L4, 100234T1 or 100234J2.

[options](#), [media options](#), [barcode](#)

From:
<https://docs.tapetrack.com/> - **TapeTrack Documentation**

Permanent link:
https://docs.tapetrack.com/options/media/relaxed_barcode_matching?rev=1573428262

Last update: **2025/01/21 22:07**

