

# Assign Volume To Slot

Volumes can be either automatically or manually assigned to a Slot.

## Automatic Slotting

Using TMSS10SlotAllocation to assign a Slot to each Volume moving to the Target Repository.

Using Restore Home Slot Repository Options.

Using Excel Import

## Manual Slotting

[Assign Slot using Volume Properties Target Repository tab.](#)

[Scan via Special Operations > Slots > Allocate to Slot](#)

[Scan Volumes to Slot Via View Repository Usage Map](#)

From:

<https://docs.tapetrack.com/> - **TapeTrack Documentation**

Permanent link:

[https://docs.tapetrack.com/technote/assign\\_to\\_slot?rev=1628648632](https://docs.tapetrack.com/technote/assign_to_slot?rev=1628648632)

Last update: **2025/01/21 22:07**

