

Lite Log Files

When completing functions in Lite, a log file is produced to record the User, Date, process and Volumes.

Log File Location

When Lite is installed, the default location for the log files is in the Documents folder of the User the program is installed under, C:\Users\user_name\Documents.

There are a couple of methods available to set the location the log file saves to if you don't want to utilize the default location.

Setting Log File via Lite


Options

When cycling through a main function in Lite (Scan-In, Scan-out, Send and Recall), after selecting the function and clicking Next you will be presented with the Options window.

To change the default directory, click in the field to the right of the Log File Folder and click on the (...) button and select the required directory to save the log file to.

This value will be recorded in the registry and held for future uses of Lite to write the log files to.

Ensure the directory you choose has the appropriate permissions for Lite to write to.



From:
<https://docs.tapetrack.com/> - **TapeTrack Documentation**

Permanent link:
https://docs.tapetrack.com/technote/log_file?rev=1599615636

Last update: **2025/01/21 22:07**

