

Labelling Promedia Drawers

Labelling the Promedia Drawers in the following manner allows the configuration of TapeTrack allowing operators to quickly locate a [Slotted Volume](#), or locate a specified [Slot](#) to place a [Volume](#) in.

With TapeTrack [Zones](#) and [Slotting](#) information configured correctly, each [Slot](#) will be referenced with the [Zone-ID](#), Drawer row and [Slot](#) number in the format of [Zone-ID](#) Row#.Slot#.

For example **Red-2 3.16** would be in the 2nd drawer (Red-2), in the third row, 16th [Slot](#).

Red-4 2.20 would be in the 4th drawer (Red-4), in the 2nd row, 20th [Slot](#) in the drawer Red-4.

From:

<https://docs.tapetrack.com/> - **TapeTrack Documentation**

Permanent link:

https://docs.tapetrack.com/technote/promedia_labels?rev=1636496676

Last update: **2025/01/21 22:07**

