

Scan To Slot Via Special Operations Slot Allocation

To use the Allocate to Slot via the Special Operations menu, the selected Volumes must be either moving to or confirmed in a Repository with allocated Zones.

While this process can be used any number of Volumes, you can only select a Slot for one Volume at a time. To Slot multiple Volumes more efficiently see [Scan to assign using Repository Usage Map](#).

Select Volumes To Slot

Select Slots

Refresh Display

From:
<https://docs.tapetrack.com/> - TapeTrack Documentation

Permanent link:
https://docs.tapetrack.com/technote/scan_to_assign_slot_allocate?rev=1628643693

Last update: **2025/01/21 22:07**

